

coreflon

PTFE SURGICAL SUTURE

PTFE SURGICAL SUTURE BLACK and WHITE

NEEDLE BLACK	needle length	needle cross- section	suture length	USP			
				3/0	4/0	5/0	6/0
DS 1/2 circle	13		55cm	-	-	C5DS135512B	C6DS135512B
	16		55cm	C3DS165512B	C4DS165512B	C5DS165512B	C6DS165512B
	19		55cm	C3DS195512B	C4DS195512B	C5DS195512B	-
	22		55cm	C3DS225512B	C4DS225512B	-	-
DKO 3/8	16		55cm	-	C4DKO165512B	C5DKO165512B	-

PTFE SURGICAL SUTURE SILVER and WHITE

NEEDLE BLACK	needle length	needle cross- section	suture length	USP			
				3/0	4/0	5/0	6/0
DS 1/2 circle	13		55cm	-	-	C5DS135512S	-
	16		55cm	C3DS165512S	C4DS165512S	C5DS165512S	-
	19		55cm	C3DS195512S	C4DS195512S	C5DS195512S	-
	22		55cm	C3DS225512S	C4DS225512S	-	-
DM 3/8 circle	13		55cm	-	-	C5DM135512S	C6DM135512S
	16		55cm	C3DM165512S	C4DM165512S	C5DM165512S	C6DM165512S
	19		55cm	C3DM195512S	C4DM195512S	C5DM195512S	-
	22		55cm	C3DM225512S	C4DM225512S	-	-

coreflon SILVER
PTFE polytetrafluoroethylene monofilament nonabsorbable surgical suture
maximum time of device contact with tissue - 30 days
product free of phthalates and APO

4/0 USP
DM16
3/8 circle
55cm

CE 0197

manufacturer
expiration date
product code
LOT number
suture material
needle colour

rev. 2
consult instruction for use
do not use if package is damaged
keep away from sunlight
do not resterilize
upper limit of temperature
do not reuse
sterilized using ethylene oxide

needle length
needle cross-section
curvature
needle type
USP size

DS: round needle with cutting blade. Provides a smooth transition between the hard tissues.

DKO: reverse cutting needle with a triangular cross section with three cutting edges where the third edge is outside curvature of the needle.

DM: needle with reverse cutting blade - the third section of the length of the needle is sharp, while the remaining part has an oval cross-section, flattened for better grip in a vise.